

S-Blade

Laryngoscope blade



Reusable Miller blade.

Specific features

Connection standard ISO 7376

Technical data

Size	4
Blade length	180 mm
Total length	205 mm
Weight	120 g
Material	Stainless steel
Light transmission	Traditional light

Class I Med.Device according to Dir. 93/42/CEE

S-BLADE SIZE 4-MILLER BLADE

DG04004A



CND Classification Z12021003

Registration number 147194

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it