

C-Blade

Macintosh blade



Macintosh blade whose bending allow a natural sliding into the position of the vallecula and the wide flange assists keeping the tongue safe when intubating.

Specific features

Connection standard ISO 7376
6 sizes available

Technical data

Sizes	0-1-2-3-4-5
Material	Stainless stel
Light	Traditional

Class I med. Dev. according to dir. 93/42/CEE

C-BLADE

DG05000A
DG05001A
DG05002A
DG05003A
DG05004A
DG05005A



CND Classification Z12021003

Rev.2 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it