

## Spare bulb



Spare bulb for traditional light laryngoscope blades.

### Specific features

Easy to replace  
Compatible with blades sizes 0-1-2-3-4

### Technical data

Overall length	19,3 mm
Diameter	5 mm
Weight	1g

TFJ SPARE LAMP SIZES 0-1-2-3-4

**DG05501B**

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)