

## ADULT EXTRA LARGE CUFF



Adult Cuff for sphygmomanometer suitable for use on the thigh.

### Specific Characteristics

Air chamber with 1 tube

### Technical data

Dimensions opened	750x220 mm
Suitable for circumferences	From 39 to 56 cm
Tube length	50 cm
Material	Nylon, rubber
Weight	215 g

THIGH ORANGE CUFF 1 TUBE W/AIR CHAMBER  
FOR SPHYGMOMANOMETER

**DG06002A**

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)