

LARGE ADULT CUFF



Adult Cuff for sphygmomanometer suitable for use on arm.

Specific Characteristics

Air chamber with 1 tube

Technical data

Dimensions opened	620x170 mm
Suitable for circumferences	From 33 to 50 cm
Tube length	50 cm
Material	Nylon, rubber
Weight	110 g

ORANGE LARGE ADULT CUFF 1 TUBE W/AIR
CHAMBER FOR SPHYGMOMANOMETER
DG06003A

Rev.2 (25/08/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it