

## PEDIATRIC CUFF



Pediatric cuff for sphygmomanometer suitable for use on arm.

### Specific Characteristics

Air chamber with 1 tube

### Technical data

Dimensions opened	340x100 mm
Suitable for circumferences	From 18,5 to 27 cm
Tube length	50 cm
Material	Nylon, rubber
Weight	45 g

PEDIATRIC ORANGE CUFF 1 TUBE W/AIR CHAMBER  
FOR SPHYGMOMANOMETER

DG06005A

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)