

Sistema 3

Multicuff sphygmomanometer



Semi-rigid case equipped with Sfigmotek, is a device for non-invasive measurement of arterial blood pressure, and 3 cuffs.

Specific Features

- Single tube measurement system
- The dial protection ensures greater impact resistance
- The large pressure gauge and the double colored graduated scale, allow easy reading
- The shape of the cuff and the metal slot, prevent from undesired opening and allow easy application
- Inflation bulb with non-slip surface
- Direction of application marked on the cuff
- The semi-rigid case offers great protection to the content
- The internal of the case has two internal separators with pockets, also marked with indication of the intended use of the stored cuff
- The kit comes with 3 cuffs: DG06003 (large adult cuff), DG06004 (Adult cuff), DG06005 (pediatric cuff)

Technical data

Case dimensions	235x285x80 ± 5mm
Case weight	630 g
Large cuff dimensions ^(opened)	600x170 mm
Large cuff circumference	From 33,3 to 51 cm
Large cuff weight	110 g
Adult cuff dimensions ^(opened)	530x140 mm
Adult cuff circumference	From 25,4 to 40,6 cm
Adult cuff weight	70 g
Pediatric cuff dim. ^(opened)	340x100 mm
Pediatric cuff circumference	From 18,5 to 27 cm
Pediatric cuff weight	45 g
Pressure gauge scale	0 – 300 mmHg
Pressure gauge diameter	73 mm
Dial diameter	60 mm
Sphygmom. weight	290 g
Materials	Al, chromed brass, Nylon

Class IM Med.Device according to Dir. 93/42/CEE

SISTEMA 3 - MULTI-CUFF
SPHYGMOMANOMETER



DG06007A

CND Classification	C9006
Registration number	882374

Rev.1 (25/08/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it