

Sistema 5

Multicuff sphygmomanometer



Semi-rigid case equipped with Sfigmotek, a device for non-invasive measurement of arterial blood pressure, and 5 cuffs.

Specific Features

- Single tube measurement system
- The dial protection ensures greater impact resistance
- The large pressure gauge and the double colored graduated scale , allow easy reading
- The shape of the cuff and the metal slot, prevent from undesired opening and allow easy application
- Inflation bulb with non-slip surface
- Direction of application marked on the cuff
- The semi-rigid case offers great protection to the content
- The internal of the case has two internal separators with pockets, also marked with indication of the intended use of the stored cuff
- The kit comes with 5 cuffs: Thigh Cuff, Large adult, Adult, Pediatric, Infant

Technical data

| | |
|---|--------------------------|
| Case dimensions | 235x285x80 ± 5mm |
| Case weight | 630 g |
| Thigh cuff dim. ^(opened) | 750x215 |
| Thigh cuff circumference | From 41 to 56 cm |
| Thigh cuff weight | 215 g |
| Large adult cuff dim. ^(opened) | 600x170 mm |
| Large adult cuff circumference | From 33,3 to 51 cm |
| Large adult cuff weight | 110 g |
| Adult cuff dim. ^(opened) | 530x140 mm |
| Adult cuff circumference | From 25,4 to 40,6 cm |
| Adult cuff weight | 70 g |
| Pediatric cuff dim. ^(opened) | 340x100 mm |
| Pediatric cuff circumference | From 18,5 to 27 cm |
| Pediatric cuff weight | 45 g |
| Infant cuff dim. ^(opened) | 250x80 mm |
| Infant cuff circumference | From 13 to 19 cm |
| infant cuff weight | 35 g |
| Pressure gauge scale | 0 – 300 mmHg |
| Pressure gauge diameter | 73 mm |
| Dial diameter | 60 mm |
| Sphygmom. weight | 290 g |
| Materials | Al, chromed brass, Nylon |

Class IM Med.Device according to Dir. 93/42/CEE

SISTEMA 5 - MULTI-CUFF
SPHYGMOMANOMETER



DG06008A

CND Classification C9006

Registration number 1530326

Rev.2 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it

Rev.2 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it