

Sistema 5

Multicuff sphygmomanometer



Semi-rigid case equipped with Sfigmotek, a device for non-invasive measurement of arterial blood pressure, and 5 cuffs.

Specific Features

- Single tube measurement system
- The dial protection ensures greater impact resistance
- The large pressure gauge and the double colored graduated scale , allow easy reading
- The shape of the cuff and the metal slot, prevent from undesired opening and allow easy application
- Inflation bulb with non-slip surface
- Direction of application marked on the cuff
- The semi-rigid case offers great protection to the content
- The internal of the case has two internal separators with pockets, also marked with indication of the intended use of the stored cuff
- The kit comes with 5 cuffs: Thigh Cuff, Large adult, Adult, Pediatric, Infant

Technical data

Case dimensions	235x285x80 ± 5mm
Case weight	630 g
Thigh cuff dim. ^(opened)	750x215
Thigh cuff circumference	From 41 to 56 cm
Thigh cuff weight	215 g
Large adult cuff dim. ^(opened)	600x170 mm
Large adult cuff circumference	From 33,3 to 51 cm
Large adult cuff weight	110 g
Adult cuff dim. ^(opened)	530x140 mm
Adult cuff circumference	From 25,4 to 40,6 cm
Adult cuff weight	70 g
Pediatric cuff dim. ^(opened)	340x100 mm
Pediatric cuff circumference	From 18,5 to 27 cm
Pediatric cuff weight	45 g
Infant cuff dim. ^(opened)	250x80 mm
Infant cuff circumference	From 13 to 19 cm
infant cuff weight	35 g
Pressure gauge scale	0 – 300 mmHg
Pressure gauge diameter	73 mm
Dial diameter	60 mm
Sphygmom. weight	290 g
Materials	Al, chromed brass, Nylon

Class IM Med.Device according to Dir. 93/42/CEE

SISTEMA 5 - MULTI-CUFF
SPHYGMOMANOMETER



DG06008A

CND Classification C9006

Registration number 1530326

Rev.1 (12/04/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it

Rev.1 (12/04/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it