

SFIGMOTEK

Palm sphygmomanometer



Sfigmotek is a device for non-invasive measurement of arterial blood pressure.

Specific Characteristics

- Single tube measurement system
- The dial protection ensures greater impact resistance
- The large pressure gauge and the double colored graduated scale , allow easy reading
- The shape of the cuff and the metal slot, prevent from undesired opening and allow easy application
- Inflation bulb with non-slip surface
- Direction of application marked on the cuff

Technical data

Cuff width	135 mm
Cuff circumference	Da 26 a 36 cm
Case dimensions	180x110x65 mm
Tube length	48 cm
Pressure gauge scale	0 – 300 mmHg
Pressure gauge diameter	73 mm
Dial diameter	60 mm
Materials	Al, chromed brass, Nylon
Weight with case	440 g

Class IM Med.Device according to Dir. 93/42/CEE

SFIGMOTEK WITH ADULT CUFF
- PALM
SPHYGMOMANOMETER



DG06022B

CND Classification C9006

Registration number 90065

Rev.1 (12/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it