

## DGX 062

### Pen Light



Light for the inspection of eye and oropharyngeal cavity.

#### Specific characteristics

Plastic body  
 Incandescent light  
 The light turns on by pressing the clip  
 Requires 2 batteries AAA not included

#### Technical data

Length	127 mm
Weight without batteries	11 g

#### Class I Med.Device according to Dir. 93/42/CEE

DGX 062 - BLUE PLASTIC PEN  
 LIGHT NO BATTERIES

DG08062D

CND Classification V9099

Registration number 785990



Rev.2 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)