

DGX 062

Pen Light



Light for the inspection of eye and oropharyngeal cavity.

Specific characteristics

Plastic body
 Incandescent light
 The light turns on by pressing the clip
 Requires 2 batteries AAA not included

Technical data

Length	127 mm
Weight without batteries	11 g

Class I Med.Device according to Dir. 93/42/CEE

DGX 062 - BLUE PLASTIC PEN
 LIGHT NO BATTERIES
 DG08062D



CND Classification	V9099
Registration number	785990

Rev.1 (09/09/2016)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it