

PEN LIGHT



Light pen to be used by professionals. Thanks to the wide beam, constitutes a valuable aid for rescuers operating in conditions of poor visibility.

Specific Features

- Requires 3 batteries LR41 included
- Led light on the upper side of the pen
- Rigid case included
- Power button on the upper end to facilitate lighting
- Ergonomic and non slipping grip

Technical data

| | |
|------------|------------|
| Dimensions | ø15x130 mm |
| Pen weight | 26 g |
| Kit weight | 41 g |

LED LIGHT PEN - LUMINOUS PEN WITH
LED
DG08080D

Rev.0 (25/09/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italy Tel. +39.0521.541111 Fax +39.0521.541222
www.spencer.it