

TEST STRIPS

For blood Glucose Meter System



Test strips for the determination of blood glucose, to be used with the instrument DGX911 - Accu-Chek.

Specific features

A larger, wider test strip is easy for everyone to use. It features a protective coating so can be handled from anywhere on the strip without risk of influencing the blood glucose test result. Accu-Chek Aviva test strips quickly attract and absorb a tiny blood drop. The wide-mouth dosing area is designed so you can fill the strip on the very first try.

Technical data

| | |
|---------------------------------------|-----------|
| Blood sample required for measurement | 0.6µl |
| Required measurement time | 5 seconds |
| Pieces included | 50 pcs |

IVD according to Dir. 93/42/CEE

REACTIVE STRIPS ACCU-CHEK AVIVA, 50 PCS-GLUCOMETER 50 PZ.



DG09114E

CND Classification

W0101060101

Distributed by Spencer Italia s.r.l.

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>

Spencer Italia S.r.l. a socio unico Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it