

DGX 212

Otoscope



DGX 212 is a medical equipment suitable to directly examine the external acoustic duct and the tympanic membrane. The Otoscope has a screw cap. DGX 212 provided with a rheostat to adjust the intensity of light, a removable glass allowing a 3X magnification and an output for air insufflation to check the mobility of the tympanic membrane.

Specific Characteristics

Reliable and indestructible high quality instrument

Invaluable for everyday ear nose and throat use

Easily adjustable Neutral 2.5V vacuum light (made in Germany)

Quality 3 times magnifying lens

Ideal for pneumatic tympanic membrane tests, can also be used for rhinoscopic examinations

Supplied with 3 autoclavable ear funnels in a colour coordinated zip bag. All latex free

Technical data

Cap	Screw type
Supplied Ear funnels	Ø 2,5 – 3,5 – 4,5 mm
Required batteries	2x C-Type not included

Class I Med.Device according to Dir. 93/42/CEE

DGX 212 - OTOSCOPE

DG10212



CND Classification

Z12149006

Distributed by Spencer Italia s.r.l.

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it