

EVX 50

Mount catheter



Flexible catheter mount for connection between breathing device and connector of tracheal tube or tracheotomy or face mask.

Specific features

Transparent tube, flexible and lightweight with smooth inner surface
 Cobb connector with tap for bronchial suction with hole for bronchoscope.
 Latex and PHT Free
 Connector 15M
 Non conductive
 The device is sterile and disposable

Technical data

Tube length	10 cm
Dead space	22 ml
Compliance	< 2% of the volume
Material	PVC
Flow resistance <small>(according to EN 12342)</small>	0,029 kpa - 15 l/min 0,078kpa - 30 l/min 0,327 kpa - 60 l/min

Medical device CE marked according to dir. 93/42/CE
 EVX 50 - MOUNT CATHETER
 DG10500K



Rev.2 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it