

## DGX 195

Disposable phonendoscope



Disposable phonendoscope for chest auscultation. The package contains 10 pieces.

### Specific Characteristics

Extra light-weight

Ideal for single-patient use in contagious areas or education purposes

### Technical data

Head diameter	48 mm
Weight	42 g/each

Class I Med.Device according to Dir. 93/42/CEE

DGX195 DISPOSABLE  
PHONENDOSCOPE 10PCS  
DG21196B



CND Classification	C9005
Registration number	90204

Rev.2 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. a socio unico Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)