

DGX 190**Sprague Rappaport stethoscope**

Sprague Rappaport stethoscope with two heads particularly versatile. The included equipment, allows to adapt the device to neonatal, pediatric and adult patients and makes it suitable for particular types of auscultations depending on the used configuration.

Specific features

Extremely light
 Two-tubes structure
 Metal clips limit the friction of the tubes
 Equipped with adult, pediatric and infant bells
 2 spare ear-tips
 Large and small spare diaphragms

Technical data

Large head diameter	45 mm
Small head diameter	31 mm
Tubes length	56 cm
Weight	235 ± 10g

Class I Med.Device according to Dir. 93/42/CEE

D - DGX 190 - SPRAGUE
 RAPPAPORT CROM/ORANGE
 DG21200C



Distributed by Spencer Italia S.r.l.

Rev.0 (29/06/2020)

COPIA NON CONTROLLATA – successive revisioni saranno disponibili al link <http://support.spencer.it>

Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it