

# DGX 30

## Intubation Stylet



Intubation stylet used to make the application of endotracheal tubes easier.

### Specific features

- Flexible aluminum core, covered by a plastic coating
- Rounded distal parts to avoid tissue damage
- The device is disposable and is provided sterile

### Technical data

Package dimensions	65x495 ± 5 mm
Weight	14 ± 4 g
Length	34 ± 1 cm
Diameter	Ø4,7 mm (14Fr)

Class IIa Med.Device according to Dir. 93/42/CEE

DGX 30 ADULT - DEFORMABLE STYLET  
DG50300A



CND Classification

L031199

Distributed by Spencer Italia S.r.l.

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)