

Pediatric cuff



Pediatric cuff for use with sphygmomanometers for non-invasive measurement of arterial blood pressure.

Specific Characteristics

System with two pipes
Cuff with indications for direction of application

Technical data

Dimensions of ext. cuff	445x100 mm
Tubes length	50 cm
Materials	Nylon, rubber
Weight	70 ± 10 g

BLUE NYLON PEDIATRIC CUFF 2
TUBES
DG60001A

Rev.1 (15/03/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it