

Pedi Intube

Intubation kit



Soft bag equipped with the devices commonly used for intubation of pediatric patients during rescue procedures.

Specific features

- Fabric made of extremely robust nylon texture
- Transparent window with strap closure on the upper side to accommodate identification plates
- Double zip closure to make opening and closure easier
- Strap band on the lower side to attach the bag inside backpacks equipped with male straps
- 5 vertical pockets with strap closure
- Internal padded separator with two open pockets and on the opposite side 6 mesh compartments
- 1 transparent pocket with zip closure

Technical data

Bag dimensions	300x110x190 ± 10 mm
Empty bag weight	270 ± 10 g
Bag volume	4,8 ± 0,2 l
Blades	Miller and Macintosh conventional light

Included equipment

DG00603A	UNIVERSAL BLACK - TUFF CUT SCISSORS
DG03101A	L-HANDLE CONVENTIONAL LIGHT PAEDIATRIC
DG04000A	S-BLADE SIZE 0-MILLER BLADE
DG04001A	S-BLADE SIZE 1-MILLER BLADE
DG05000A	C-BLADE - SIZE 0-MACINTOSH BLADE
DG05001A	LARYNGOSCOPE BLADE TYPE MCINTOSH SIZE 1
DG80531K	C-F UN-CUFFED TRACHEAL TUBE SIZE 3
DG80536K	C-F UN-CUFFED TRACHEAL TUBE SIZE 3,5
DG80541K	C-F UN-CUFFED TRACHEAL TUBE SIZE 4
FG00607A	FORCEPS MAGILL PEDIATRIC
GM03050K	STERILE DISPOSABLE SYRINGE 5ML
GM03100K	DISPOSABLE STERILE SYRINGE 10
GM03225K	SILK BANDAID 2,5cm x 9,14m
RM10320A	GUEDEL SIZE 50 MM
RM10321A	GUEDEL SIZE 60 MM
RM10325A	GUEDEL SIZE 40 MM
RM30700A	MOUNDER - JAW SPREADER

Med. Dev. CE marked according to dir. 93/42/CE

PEDI INTUBE COMPLETE - PAEDIATRIC
SOFT INTUBATION
DG80007B



Rev.2 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it