

MATRIX INTUBE

Professional intubation Kit



Semi-rigid bag equipped with the devices commonly used for intubation during rescue operations.

Included equipment	
DG00603A	UNIVERSAL BLACK - TUFF CUT SCISSORS
DG03100A	L-HANDLE CONVENTIONAL LIGHT ADULT
DG04001A	S-BLADE - MILLER BLADE SIZE 1
DG04002A	S-BLADE - MILLER BLADE SIZE 2
DG04003A	S-BLADE - MILLER BLADE SIZE 3
DG05001A	C-BLADE - MACINTOSH BLADE SIZE 1
DG05002A	C-BLADE - MACINTOSH BLADE SIZE 2
DG05003A	C-BLADE - MACINTOSH BLADE SIZE 3
DG50300C	DGX 30 ADULT - DEFORMABLE STYLET
DG80550K	C-F CUFFED TRACHEAL TUBE SIZE 5
DG80560K	C-F CUFFED TRACHEAL TUBE SIZE 6
DG80570K	C-F CUFFED TRACHEAL TUBE SIZE 7
DG80580K	C-F CUFFED TRACHEAL TUBE SIZE 8
DG90550K	DGX 550 - WENDL NASOPHARYNG. CANN.SZ.0 (5 mm)
DG90560K	DGX 550 - WENDL NASOPHARYNG. CANN.SZ.2 (6 mm)
DG90570K	DGX 550 - WENDL NASOPHARYNG. CANN.SZ.4 (7 mm)
FG00605A	MAGILL ADULT - TUBE-GUIDING FORCEPS
FG00607A	FORCEPS MAGILL PEDIATRIC
GM03050K	GMX05 - STERILE DISPOSABLE SYRINGE 5ML
GM03100K	GMX 10 - STERILE DISPOSABLE SYRINGE 10 ml
GM03225K	S-PLASTER - SILK BANDAID 2,5cm x 9,14m
RM10320A	GUEDEL SIZE 50 mm
RM10321A	GUEDEL SIZE 60 mm
RM10322A	GUEDEL SIZE 70 mm
RM10323A	GUEDEL SIZE 80 mm
RM10324A	GUEDEL SIZE 90 mm
RM10325A	GUEDEL SIZE 40 mm
RM30700A	MOUNDER - JAW SPREADER

Specific features	
Semi-rigid case with zip closure and carrying handle	
Reflective inserts to increase visibility	
Two semi-rigid separators for a better organization of the content	
Elastic bands on detachable supports anchored to the separators by means of strap inserts	
Traditional light laryngoscope	

Technical data	
Bag dimensions	310x360x100±5mm
Empty bag weight	850 ± 10 g
Materiali	EVA
Detachable supports	8 of different sizes

Med. Dev. CE marked accordint to dir. 93/42/CE	
MATRIX INTUBE COMPLETE - SEMI-RIGID	
PROFESSIONAL INTUBATION KIT	
DG80008B	



Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it