

# RAKAIA INTUBE

## Intubation kit



Leg waist bag ideal to face hardest work conditions, where the use of more bulky products would be difficult. In this version is equipped with devices most commonly used for intubation during rescue procedures.

### Specific features

- Extremely resistant fabric
- Upper belt with strap to attach the pouch on the belt
- 2 belts for the leg with quick release buckle
- Frontal pocket with strap closure and quick release buckle with 3 separate compartments, one of which with elastic bands for fastening of accessories.
- Adjacent to the front pocket there is a compartment with access from both sides, closed by means of two flaps with strap closure
- Large pocket with access from the top with additional pocket without closure

### Included equipment

DG00603A	UNIVERSAL BLACK - TUFF CUT SCISSORS
DG03100A	L-HANDLE CONVENTIONAL LIGHT ADULT
DG04000A	S-BLADE - MILLER BLADE SIZE 0
DG04001A	S-BLADE - MILLER BLADE SIZE 1
DG04002A	S-BLADE - MILLER BLADE SIZE 2
DG04003A	S-BLADE - MILLER BLADE SIZE 3
DG04004A	S-BLADE - MILLER BLADE SIZE 4
DG05000A	C-BLADE - MACINTOSH BLADE SIZE 0
DG05001A	C-BLADE - MACINTOSH BLADE SIZE 1
DG05002A	C-BLADE - MACINTOSH BLADE SIZE 2
DG05003A	C-BLADE - MACINTOSH BLADE SIZE 3
DG05004A	C-BLADE - MACINTOSH BLADE SIZE 4
DG80531K	C-F UN-CUFFED TRACHEAL TUBE SIZE 3
DG80536K	C-F UN-CUFFED TRACHEAL TUBE SIZE 3,5
DG80541K	C-F UN-CUFFED TRACHEAL TUBE SIZE 4
DG80545K	C-F - CUFFED TRACHEAL TUBE SIZE 4,5
DG80550K	C-F - CUFFED TRACHEAL TUBE SIZE 5

DG80560K	C-F - CUFFED TRACHEAL TUBE SIZE 6
DG80570K	C-F - CUFFED TRACHEAL TUBE SIZE 7
DG80580K	C-F - CUFFED TRACHEAL TUBE SIZE 8
FG00605A	MAGILL ADULT - TUBE-GUIDING FORCEPS
FG00607A	FORCEPS MAGILL PEDIATRIC
GM03050K	GMX05 - STERILE DISPOSABLE SYRINGE 5ML
GM03100K	GMX 10 - STERILE DISPOSABLE SYRINGE 10 ml
GM03225K	S S-PLASTER - SILK BANDAID 2,5cm x 9,14m
RM10320A	GUEDEL - DISPOSABLE CANNULA 50 mm
RM10321A	GUEDEL - DISPOSABLE CANNULA 60 mm
RM10322A	GUEDEL - DISPOSABLE CANNULA 70 mm
RM10323A	GUEDEL - DISPOSABLE CANNULA 80 mm
RM10324A	GUEDEL - DISPOSABLE CANNULA 90 mm
RM10325A	GUEDEL - DISPOSABLE CANNULA 40 mm
RM30700A	MOUNDER - JAW SPREADER

### Technical data

Dimensions	300x220 ±10mm
Weight <sup>(empty waist bag)</sup>	300 ±10 g
Leg belts adjustable	From 30 to 100 cm
Material	Nylon

Med. Dev. CE marked accordint to dir. 93/42/CE

RAKAIA INTUBE COMPLETE - COMPACT INTUBATION KIT WITH LEG FIXATION

DG80009B

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)