

PRO INTUBE

Professional intubation kit



Semi-rigid bag equipped with traditional light laryngoscope, cannula and accessories used during intubation procedures.

Specific features

Robust nylon fabric
Transparent window with strap on the upper side to apply identification labels
Closure with double zip to make closure and opening easier
Strap on lower side to allow fixation of the bag
2 transversal transparent pockets with zip closure
1 transparent pocket with zip closure

Technical data

Dimensions	300x110x190 ± 10 mm
Bag weight	260 ± 10 g
Material	Nylon

Equipment

DG03100A	ADULT LARINGOSCOPE TRADITIONAL LIGHT
DG03103B	GRIPPER ERGONOMIC HANDLES
DG05001A	C-BLADE – MACINTOSH BLADE SIZE.1
DG05002A	C-BLADE - MACINTOSH BLADE SIZE.2
DG05003A	C-BLADE - MACINTOSH BLADE SIZE.3
DG80550K	C-F- CUFFED TRACHEAL TUBE SZ 5
DG80560K	C-F- CUFFED TRACHEAL TUBE SZ 6
DG80570K	C-F- CUFFED TRACHEAL TUBE SZ 7
DG80580K	C-F- CUFFED TRACHEAL TUBE SZ 8
FG00605A	MAGILL – ADULT TUBE-GUIDING FORCEPS 268 mm
FG00607A	MAGILL – PEDIATRIC TUBE GUIDING FORCEPS 200 mm
GM03215K	S PLASTER TAPE ADHESIVE PLASTER 1,25cmx9,14m
GM03050K	GMX05 – STERILE DISPOSABLE SYRINGE 5ml
RM30700A	MOUNDER MOUTH OPENER
RM10300A	GUEDEL – DISPOSABLE CANNULA 50 mm
RM10301A	GUEDEL - DISPOSABLE CANNULA 60 mm
RM10302A	GUEDEL - DISPOSABLE CANNULA 70 mm
RM10303A	GUEDEL - DISPOSABLE CANNULA 80 mm
RM10304A	GUEDEL - DISPOSABLE CANNULA 90 mm

PRO INTUBE – PROFESSIONAL
INTUBATION KIT W/TRADITIONAL LIGHT
BLADES
DG80010A

Rev.0 (04/09/2018)

COPIA NON CONTROLLATA – successive revisioni saranno disponibili al link <http://support.spencer.it>

Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it