

## C-F

### Cuffed tracheal tube



Device for tracheal intubation and ventilation equipped with inflatable cuff.

#### Specific features

Connection ISO 15mm

The use of thermosensitive material increases flexibility when comes in contact with body temperature

Body with radiopaque markings related to size and depth of insertion

The cuff allows perfect adhesion to the tracheal wall preventing the entry of fluids in the lungs

Single sterile package (valid 5 years from production)

The device is disposable and not re-sterilizable

Contains DEHP phtalates

#### Technical data

Material	PVC
Sterilization method	EO
Size (mm)	5

Class IIa Med.Device according to Dir. 93/42/CEE

C-F - ENDOTRACHEAL STERILE

TUBE SZ 5

DG80550K



CND Classification

R01030201

Distributed by Spencer Italia s.r.l.

Rev.1 (16/04/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[spencer.it](http://spencer.it)