

C-F

Cuffed tracheal tube kit



Device for tracheal intubation and ventilation equipped with inflatable cuff. The kit includes 10 sizes.

Specific Characteristics

Connection ISO 15mm

The use of thermosensitive material increases flexibility when comes in contact with body temperature

Body with radiopaque markings related to size and depth of insertion

The cuff allows perfect adhesion to the tracheal wall preventing the entry of fluids in the lungs

Single sterile package (valid 5 years from production)

The device is disposable and not re-sterilizable

Contains DEHP phthalates

Technical data

Material	PVC
Sterilization method	EO
Included sizes (mm)	4,5/5/6/5,5/6,5/7/7,5/8/8,5/9

Class IIa Med.Device according to Dir. 93/42/CEE

C-F - ENDOTRACHEAL STERILE

TUBE KIT 10 SIZES

DG80591K



CND Classification

R01030201

Distributed by Spencer Italia s.r.l.

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it