

DGX 550

Nasopharyngeal airway



Disposable nasopharyngeal airway designed to be inserted into nasal passageway to secure an open airway.

Specific Characteristics

Provided sterile

Technical data

Material	PVC
Internal diameter	5,5 mm
Included airways	10 pcs

Class IIa Med.Device according to Dir. 93/42/CEE

DGX 550 SIZE 1 - 5,5 MM, 10 PCS-
NASOPHARYNGEAL CAN

DG90557K



Distributed by Spencer Italia s.r.l.

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it