

DGX 550

Nasopharyngeal airway



Disposable nasopharyngeal airway designed to be inserted into nasal passageway to secure an open airway.

Specific Characteristics

Provided sterile

Technical data

Material	PVC
Internal diameter	7 mm
Included airways	10 pcs

Class IIa Med.Device according to Dir. 93/42/CEE

DGX 550 SIZE 4 - 7,0 MM, 10 PCS-
NASOPHARYNGEAL CAN

DG90571K



Distributed by Spencer Italia s.r.l.

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it