

## N-GLOVEX

Nitril gloves



N-Gloves are protective gloves, ideal to be used for examinations or first in first air procedures.

### Specific features

Do not contain latex and powder  
 Disposable and not sterile  
 Ambidextrous  
 They adapt perfectly to the shape of the hand

### Accessories

**FG04574** - RDG - RDG YELLOW - GLOVES DISPENSER

### Compliant with

ASTM D6319/00

### Technical data

Package dimensions	240x120x60 mm
Overall weight	480 ± 20g
Size	L
N° of pieces	100 pz

Class I Med.Device according to Dir. 93/42/CEE

N-GLOVEX SIZE L - NITRIL GLOVES



**DR10372K**

**CND Classification** T010299

Rev.2 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l a socio unico. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[spencer.it](http://spencer.it)