


Dolphin

Universal linear life buoy



Linear life buoy to be used as buoyancy aid of victim during water rescue. Placed around the victim's abdomen, allows the floating in supine position.

Specific features

- The segmented body confers easy grip and allow an easy placing around the victim's body
- The tow strap integrated in the main body, ensures high reliability
- The width of the shoulder strap tape, provide more comfort to the rescuer
- Metal carabiner and closure rings
- The handles offer great versatility
- The built-in face mask allows to carry out primary resuscitation procedures

Standard equipment

TA09050A flex mask w/box blue logo spencer

Optional

DY32401B HIGH VISIBILITY STORAGE BAG WITH REFLECTIVE TAPE

Technical data

Main body dimensions	1200x175x85 ± 10 mm
Tow strap length	2200 mm
Belts width	30mm / 50mm shoulder strap
Materials	PE esp, Nylon
Weight	2,1 ± 0,1 kg
Minimum buoyancy force	100 N

Class I med. Dev. According to Dir. 93/42/CEE

DOLPHIN - SALVAGENTE
LIN.UNIV.C/FLEX M.E MANIGLIE
DY32500D



CND Classification V9099
Registration number 887374

Rev.1 (03/02/2015)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>

Spencer Italia S.r.l. a socio unico Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it