

ENDURO

Multilevel stretcher



Spencer ENDURO is a multilevel stretcher for transport of patients characterized by an X-Frame structure which makes the device extremely robust. It can be used for rescue operations on road and sanitary vehicles.

Specific features

- Polyethylene board easy to clean
- Adjustable backrest
- Trendelenburg antishock position
- Adjustable heights system capable to reach 9 different positions
- Four pivoting wheels to increase sliding action
- Brakes on all 4 wheels
- Yellow painted frame
- Safety system for stretcher unloading

Accessories

| | |
|----------|--|
| EN90003C | END-T Table for stretchers 10g cert. |
| ST00491A | STX 90 Telescopic head rest for stretchers |
| ST42100A | FP – REAR FASTENER |
| ST42200A | FA – FRONT FASTENER |
| CB09025C | Tanker oxygen tank holder for stretcher side bar |
| EN90019A | AMERICA 38 FIXING SYSTEM (REQUIRES EN90000) |
| EN90000B | AMERICA 34 FASTNER FOR ENDURO |
| EN90001 | SAFETY SYSTEM FOR UNLOADING |
| IF01055B | TRACK 5 telescopic iv pole for Enduro |

Standard equipment

| | |
|----------|---|
| ST70002A | STX 702 two pieces black reflex metal belts |
|----------|---|

Technical Data

| | | |
|------------------|--|-----------------|
| Length | 1970 mm | |
| Shortened length | 1560/1660 mm | |
| Width | 570 mm | |
| | Head end | Feet end |
| Height 1 | 1060 mm | 830 mm |
| Height 2 | 990 mm | 820 mm |
| Height 3 | 910 mm | 800 mm |
| Height 4 | 830 mm | 770 mm |
| Height 5 | 750 mm | 740 mm |
| Height 6 | 660 mm | 690 mm |
| Height 7 | 570 mm | 640 mm |
| Height 8 | 460 mm | 550 mm |
| Height 9 | 340 mm | 450 mm |
| | <i>Measurement referred to the patient board</i> | |
| Weight | 40 kg | |
| Loading capacity | 250 kg | |
| Wheels | ø 200 mm | |

Class I Med.Device according to Dir. 93/42/CEE

ENDURO BLACK BOARD YELLOW FRAM
E MULTI LEVEL STRET
EN10054



CND Classification

V0899

Rev.1 (06/03/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it