

ENDURO

Multilevel stretcher



Spencer ENDURO is a multilevel stretcher for transport of patients characterized by an X-Frame structure which makes the device extremely robust. It can be used for rescue operations on road and sanitary vehicles.

Specific features

Polyethylene board easy to clean
Adjustable backrest
Trendelenburg antishock position
Adjustable heights system capable to reach 9 different positions
Four pivoting wheels to increase sliding action
Brakes on all 4 wheels
Yellow painted stainless steel frame
Safety system for stretcher unloading

Accessories

EN90003C	END-T Table for stretchers 10g cert.
ST00491A	STX 90 Telescopic head rest for stretchers
ST42100A	FP – REAR FASTENER
ST42200A	FA – FRONT FASTENER
CB09025C	Tanker oxygen tank holder for stretcher side bar
EN90019A	AMERICA 38 FIXING SYSTEM (REQUIRES EN90000)
EN90000B	AMERICA 34 FASTNER FOR ENDURO
EN90001	SAFETY SYSTEM FOR UNLOADING
IF01055B	TRACK 5 telescopic iv pole for Enduro

Standard equipment

ST70002A	STX 702 two pieces black reflex metal belts
----------	---

Technical Data

Length	1970 mm	
Shortened length	1560/1660 mm	
Width	570 mm	
	Head end	Feet end
Height 1	1060 mm	830 mm
Height 2	990 mm	820 mm
Height 3	910 mm	800 mm
Height 4	830 mm	770 mm
Height 5	750 mm	740 mm
Height 6	660 mm	690 mm
Height 7	570 mm	640 mm
Height 8	460 mm	550 mm
Height 9	340 mm	450 mm
	<i>Measurement referred to the patient board</i>	
Weight	40 kg	
Loading capacity	300 kg	
Wheels	ø 200 mm	

Class I Med.Device according to Dir. 93/42/CEE

ENDURO 300 BLACK BOARD YELLOW
FRAM E MULTI LEVEL STRET
EN10055



CND Classification

V0899

Rev.1 (29/06/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it