


AMERICA 34

Fastener for Enduro



Locking pin to be mounted on the lower frame of Enduro, It allows to fix stretcher Enduro into the ambulance when equipped with fastening system America 35 (EN90010B).

Technical data	
Dimension	120x40x65±5 mm
Material	Aluminum / Steel
Weight	1 kg approx

Class I Med.Device according to Dir. 93/42/CEE	
AMERICA 34 – FASTENING FOR ENDURO	
	
EN90000A	
CND classification	V9099

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222
www.spencer.it