

END-T

Instruments table for stretchers



Table for instruments and tools to be mounted on the frame of Spencer stretchers Cross, Carrera and Enduro. The table can be used to support monitors, defibrillators and accessories.

Specific Characteristics

- Compact and lightweight
- Customizable thanks to the MARK UP system
- The structure allows quick installation and ensures easy and safe 10g fastening
- Strap bands allow easy positioning of optional bags or instruments
- Belt with integrated reel

Accessories

CB09025C	Tanker - oxygen tank holder
EN90012C	IV pole for END-T
ST02106A	Mark up – Perm. Personalisation system

Compliant with

EN 1789

Technical data

Useful support surface	375x590 mm
Thickness folded (without accessories)	90 mm
Distance between stretcher's patient floor and the bottom part of the table ¹	295 mm
Dimensions folded	400x620x90 mm
Maximum belt length (useful)	820 mm
Loading capacity referred to EN1789	15 kg
Loading capacity not referred to EN1789	25kg
Frame material	Al
Support surface material	PE
Weight	4,6 kg

¹ Can be different depending on the stretcher used.

Class I Med.Device according to Dir. 93/42/CEE

END-T TRAY FOR STRETCHER 10G
CERTIFIED
EN90003C



CND Classification	V9099
Registration number	416482

Rev.5 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it