


AMERICA 38

Fastener for Enduro



Fastening system for Enduro stretcher included a security system for the unloading device. The device is to be permanently installed inside the ambulance.

Specific features

- Reduced dimension
- Simple and intuitive fastening system
- The security system prevents an uncontrolled descent of the stretcher by ambulance

Technical data

Front fastening dimension	430x455x75 mm
Lateral fastening dimension ^(Close)	780x120x260
Material	Steel
Weight	4 kg

Class I Med.Device according to Dir. 93/42/CEE

AMERICA 38 – FIXING SYSTEM FOR ENDURO



EN90019A

CND classification

V9099

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222
www.spencer.it