


PATIENT CIRCUIT



Disposable patient circuit intended for use with Spencer lung ventilators.

Specific features

Tube connection 22F
Inner side of the tube is smooth to avoid stagnation, outer side with spirals to guarantee an inner constant lumen
High bacterial and viral protection
Filter with Luer lock port
Integrated overpressure valve
Disposable
Latex Free

Technical data

Materials	PVC, PC
Weight	200 ± 10g
BEF and VEF efficiency	>99,99%
Overpressure valve	60 cmH ₂ O
Tube length	1,2 m

Ass. Med.dev. according to dir. 93/42/CEE

PATIENT CIRCUIT
EV00106A



Registration number 4991

Rev.2 (18/08/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it