

O-Tube

Oxygen connection tube



Tube intended for pneumatic connection of devices requiring the use of pressurized oxygen.

Specific features

Colored and marked for use with oxygen
 Ideal for low pressure applications (max 14bar)
 DEHP FREE

Technical data

Internal diameter	6,7 mm
External diameter	12,7 mm
Length	100 cm
Material	fabric reinforced PVC
Burst pressure	56 bar/23°C
Output connection	Quick coupling socket
Input connection	Threaded M1/4" G

DRIVING OXYGEN HOSE CM.
 100
EV30020A

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>

Spencer Italia S.r.l. a socio unico Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it