

EVX 07

Non-rebreathing valve



Non-rebreathing valve ideal for use with manual resuscitators.

Specific features

Connection standard ISO 15/22 mm
Silicone membranes
Autoclavable

Technical data

Dimensions	87x50x85 mm
Overpressure valve	40cmH ₂ O
Materials	PC, Silicone
Weight	45 g

Class IIa med. dev. according to Dir. 93/42/CEE

EVX 07 - OVERPR.VALVE
NON/REB.AUTOCLAV.PAED./
NEON.

EV50007C



CND Classification	R0280
Registration number	310100

Rev.2 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it