

EVX 12

Non-rebreathing valve



Non-rebreathing valve ideal for use with manual resuscitators.

Specific features

Connection standard ISO 15/22 mm
Silicone membranes
Autoclavable

Technical data

Materials	PC, Silicone
Weight	40 g

Class IIa med. dev. according to Dir. 93/42/CEE

EVX 12 AUTOCLAVABLE -
NON-REBREATHING VALVE
FOR MAN. RESUSCITATOR
EV50012C



CND Classification	R0280
--------------------	-------

Rev.0 (05/09/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it