

10G WALL BRACKET



10G Wall bracket for Kompak series lung ventilators..

Specific Characteristics

Polished surfaces
Fastener with quick coupling system
Easy installation thanks to the holes that allow easy leveling

Compliant with

UNI EN 1789

Technical data

Materials	Stainless steel, Nylon
Dimensions	415x465x165 mm
Weight	5,7 kg

Class I Med.Device according to Dir. 93/42/CEE

10G WALL FIXATION SYSTEM
FOR KOMPAK

EV50107E



CND Classification V9099

Registration number 1247645

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it