

NXT MAX

10 G Fastener for lung ventilators NXT



NXT MAX fixing systems are intended for installation in ambulances to anchor Spencer ventilators.

Specific features

Its shape is designed to allow easy and quick placing of the lung ventilator on the fastener

The quick release button on the front side is easy to operate

Thanks to the perfect integration between the fastener and the lung ventilator, the overall dimensions remain almost the same as those of the lung ventilator alone, increasing its height by only a couple of centimeters

Easy to install inside the emergency vehicle, is fastened using the base

Compatible with 118 NXT, 170 NXT, 190 NXT lung ventilators

Technical data

Dimensions (DxLxH)	160 x 160 x 95 ± 5 mm
Materials	HDPE, Nylon, Steel
Weight	350 ± 10 g

Reference standard

EN 1789:2007 + A2:2014

Class I MD compliant with UE Reg. 2017/745

NXT MAX - 10G QUICK RELEASE
FASTENER
EV50120B



Rev.0 (30/11/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza(PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it