

NXT MAX WALL

10 G Wall fastener for NXT lung ventilators



NXT MAX Wall fixing systems are intended for installation in ambulances to anchor Spencer ventilators.

Specific features

Its shape is designed to allow easy and quick placing of the lung ventilator on the fastener

Thanks to the perfect integration between the fastener and the lung ventilator, the overall dimensions remain almost the same as those of the lung ventilator alone, increasing its height by only a couple of centimeters

Easy to install inside the emergency vehicle, is fastened using the back side, ideal solution for space optimization

Compatible with 118 NXT, 170 NXT, 190 NXT lung ventilators

Technical data

Dimensions (DxLxH)	160 x 160 x 95 ± 5 mm
Materials	HDPE, Nylon, Steel
Weight	300 ± 10 g

Reference standard

EN 1789:2007 + A2:2014

Class I MD compliant with UE Reg. 2017/745

NXT MAX WALL - 10G WALL
FASTENER
EV50122B



Rev.0 (30/11/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza(PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it