

MAGILL

Tube-guiding forceps



Pinza of Magill useful during the phases of an adult patient intubation.

Specific features

- Sterilizable
- Reeded distal part to offer a better grip.

Technical data

Total length	245 mm
Curved part length ¹	135 mm
Weight	77 g
Material	Stainless steel

¹ From distal part to first curving

Class I Med.Device according to Dir. 93/42/CEE

FORCEPS MAGILL ADULT

FG00605A



CND Classification L031312

Registration number 141543

Rev.2 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it