

MAGILL

Tube-guiding forceps



Magill forceps useful during intubation procedures of pediatric patients.

Specific features

Sterilizable
 Reeded distal part to offer a better grip.
 Satin finish

Technical data

Total length	200 mm
Curved part length ¹	100 mm
Weight	46 g
Material	Stainless steel

¹ From distal part to first curving

Class I Med.Device according to Dir. 93/42/CEE

FORCEPS MAGILL CHILDREN
 FG00607A



CND Classification	L031312
Registration number	141547

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it