

# GMX 295

## Bandage scissors



Scissors to be used for cut of bandages or clothes.

### Specific features

- Sterilizable
- The build quality allows a precise cut
- The shape of the end allows easy insertion under clothes or bandages.

### Technical data

Overall length	184 mm
Blade length	45 mm
Weight	78 g
Material	Stainless steel

Class I Med.Device according to Dir. 93/42/CEE

LISTER BANDAGE SCISSORS

FG02295A



CND Classification L031305

Registration number 135992

Rev.2 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[spencer.it](http://spencer.it)