

FGX 320
Buck hammer



Reflex hammer to test deep tendon reflexes.

Specific characteristics

- 2 Heads with different dimensions
- Equipped with integrated needle and brush

Technical data

Length	180 mm
Heads dimensions	ø14 and 12 mm
Weight without batteries	85 g

FGX 320 - BUCK HAMMER



FG02320B

Rev.2 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it