

GMX 308

Surgical forceps



Stainless steel surgical forceps.

Specific features

- Sterilizable
- Reeded I part to offer a better grip.
- The absence of screws, makes it very stable and reliable

Technical data

Total length	140 mm
Material	Stainless steel

Class I Med.Device according to Dir. 93/42/CEE

GMX 308 - SURGICAL
FORCEPS

FG03308B

CND Classification L031399

Registration number 135987



Rev.1 (07/06/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it