

# GMX 350

## Kelly forceps



Kelly curved hemostatic forceps.

### Specific features

- Sterilizable
- Reeded distal part to offer a better grip.
- The absence of screws, makes it very stable and reliable

### Technical data

Total length	145 mm
Clamping area length	30 mm
Reeded area length	20 mm
Weight	32 g
Material	Stainless steel

Class I Med.Device according to Dir. 93/42/CEE

KELLY CURVED HEMOSTATIC  
FORCEPS

**FG03350B**

**CND Classification**



L031306

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>

Spencer Italia S.r.l a socio unico. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[spencer.it](http://spencer.it)