

GMX 351

Kelly forceps



Kelly straight hemostatic forceps..

Specific features

- Sterilizable
- Reeded distal part to offer a better grip.
- The absence of screws, makes it very stable and reliable

Technical data

Total length	143 mm
Clamping area length	30 mm
Reeded area length	19 mm
Weight	27 g
Material	Stainless steel

Class I Med.Device according to Dir. 93/42/CEE

KELLY STRAIGHT
HEMOSTATIC FORCEPS
FG03351B



CND Classification	L031306
Registration number	128728

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l a socio unico. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it