

## RDG

### Gloves dispenser



Dispenser for table or wall mounting to store a gloves box up to to 100 pieces.

#### Specific features

- The smooth surfaces allow easy cleaning
- Rounded edges to ensure operator safety
- Easy removable cover for replacement of the gloves box contained inside
- Space saver

#### Technical data

Dimensions	265x140x95 mm
Material	PC
Weight	310 g

RDG TRANSPARENT - GLOVES  
DISPENSER  
**FG04575C**

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>

Spencer Italia S.r.l a socio unico. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[spencer.it](http://spencer.it)